

ZĒNO

It's revolutionary.
It's fast. It's clinically
proven to get rid of pimples.



ACTUAL SIZE

ZĒNO

For more information about Zeno,
call us or visit myzeno.com.

Tyrell, Inc.
Houston, Texas 77067
888-4MY-ZENO or 888-469-9366

PRINTED IN THE U.S.A. 02/05

The names Tyrell, Zeno and ClearPoint, the
Tyrell and Zeno logos, and the phrase Clearly
Outsmarts Pimples, are all trademarks,
service marks, registered trademarks or trade
names of Tyrell, Inc. Zeno is covered
by one or more of the following patents:
U.S. Patents 6,635,075; 6,245,093.



ZĒNO™

Clearly outsmarts pimples.™

In an FDA-reviewed, controlled clinical trial, 90% of blemishes treated with Zeno disappeared or faded within 24 hours.



BEFORE

AFTER

What is Zeno?

Zeno is a handheld, portable electronic medical device that is clinically proven to make pimples disappear fast. In fact, for treating acne pimples, it's the most scientifically advanced and effective device available without a prescription.

What does it do?

Zeno applies a precisely controlled heat dose directly to the pimple through a metal pad. One treatment lasts 2½ minutes. Two to three treatments spread over 24 hours clear up most pimples quickly. It's that simple.

How does it work?

Pimples occur when excess skin oil and flaking skin plug the opening of a hair follicle and create the ideal environment for acne bacteria to grow. Precisely controlled heat applied for a specific period of time causes the bacteria to self-destruct. When the bacteria die, the pimple goes away.

How many treatments does it take?

Two to three 2½-minute treatments spread over 24 hours are sufficient for most pimples. Often, only one treatment is needed.

Who should use Zeno?

Zeno is for people with mild to moderate acne. If in doubt, ask your doctor.

When should Zeno be used?

Use Zeno at the first sign of a pimple. It works best when used on a pimple in the early stages of development. Also, the sooner pimples are cleared up, the less chance they have to leave scars.

How does Zeno compare with topical or oral treatments?

Zeno can clear up a pimple in just hours. Most other medications can take days or even weeks to work. What's more, other than slight, temporary skin redness, Zeno has none of the side effects described by most topical or oral medications.

Is Zeno safe and painless?

Patients included in our clinical trials experienced little reaction to the low-level heat. Other than some temporary skin redness following a treatment cycle, no other adverse side effects were reported.

Is Zeno all I need for acne?

Zeno does not replace good skin hygiene or a doctor's care. It is designed to treat mild to moderate acne. Your doctor may recommend other therapies, depending on the severity of your acne. Follow your doctor's advice and use Zeno for the occasional pimple.

How do I learn more about Zeno?

Ask your doctor whether Zeno is right for you or visit our web site at myzeno.com.